


Introduction

Biochip array technology enables multi-analytical screening of drugs of abuse from a single undivided sample in under 20 minutes with the new biochip analyser Evidence MultiSTAT. This study reports the analytical evaluation of this application to the simultaneous screening of drugs of abuse in urine (creatinine included) and oral fluid.

Methodology

The DOA Array Urine and DoA Array Oral Fluid II were used (EV4193, EV4279, Randox Toxicology Ltd., Crumlin, UK), simultaneous competitive chemiluminescent immunoassays on a biochip surface were applied to the fully automated Evidence MultiSTAT analyser, which processes a self-contained cartridge containing all the components required for the assays (EV4115, Randox Toxicology Ltd, Crumlin, UK). Sampling against a cut-off sample, the results are qualitative. The NeoSal Oral Fluid Collection System was used for oral fluid (Neogen Corporation, Lansing, USA).



Test Menu and Cut-Offs

DOA Array Urine

Analyte	Cut-offs (Urine)	
	Cut-Off (ng/mL)	Analyte
AB-PINACA	2.5	JWH-018
α-PVP	5	6-MAM
Amphetamine	200	Metadone
Barbiturates	200	Methamphetamine
Benzodiazepines I	150	Opiate
Benzodiazepines II	150	Oxycodone
Benzoyllecgonine/Cocaine	150	THC
Buprenorphine	1	Tramadol
ETG	750	Tricyclic Antidepressants (TCA)
Fentanyl	2	UR-144
Creatinine	20 mg/dL	

DOA Array Oral Fluid II

Analyte	Cut-offs (Oral Fluid)	
	Cut-Off (ng/mL)	Analyte
α-PVP	2.5	LSD
Amphetamine	60	6-MAM
Barbiturates	60	Metadone
Benzodiazepines I	15	Methamphetamine
Benzodiazepines II	15	Opiate
Benzoyllecgonine/Cocaine	30	Oxycodone
Buprenorphine	1.5	PCP
Fentanyl	1.5	Tramadol
JWH-018	20	THC
Ketamine	65	UR-144

Results

DOA Array Urine

Repeatability

Repeatability was determined by assessing control material prepared at the cut-off and at ±50% of the cut-off. Each sample was assessed against the cut-off material twice a day for 10 days, resulting in n=20 results for each sample. The percentage agreement was calculated for the number of samples that correctly reported negative and positive.

Assay	-50% cut-off	Cut-off	+50% cut-off	Agreement (%)
AB-PINACA	0	8	19	97.5
α-PVP	0	10	20	100
Amphetamine	0	7	20	100
Barbiturates	0	8	20	100
Benzodiazepines I	0	12	20	100
Benzodiazepines II	0	11	20	100
Benzoyllecgonine/Cocaine	0	12	20	100
Buprenorphine	0	8	20	100
ETG	0	12	1	97.5
Fentanyl	0	13	20	100
JWH-018	0	12	20	100
6-MAM	0	14	20	100
Metadone	0	12	20	100
Methamphetamine	0	12	20	100
Opiate	0	12	20	100
Oxycodone	0	7	20	100
THC	0	14	20	100
Tramadol	0	12	20	100
TCA	0	11	20	100
UR-144	0	9	20	100
Creatinine	0	10	20	100

Accuracy

Accuracy was determined by assessing spiked samples at varying concentrations (50 spiked positive samples prepared at concentrations greater than the cut-off, 10 negative spiked samples prepared at concentrations lower than the cut-off and 40 blank negative samples). Each sample was assessed against the cut-off material to determine a positive or negative result. The percentage agreement was calculated as the percentage of correct reports out of the total number of samples analysed (n=100).

Assay	Spike +	Spike -	Agreement (%)
AB-PINACA	50	0	100
α-PVP	50	0	100
Amphetamine	50	0	100
Barbiturates	50	0	100
Benzodiazepines I	50	0	100
Benzodiazepines II	50	0	100
Benzoyllecgonine/Cocaine	50	0	100
Buprenorphine	46	0	96
ETG	50	0	100
Fentanyl	50	0	100
JWH-018	50	0	100
6-MAM	50	0	100
Metadone	50	0	100
Methamphetamine	50	0	100
Opiate	50	0	100
Oxycodone	50	0	100
THC	44	0	94
Tramadol	50	0	100
TCA	50	0	100
UR-144	50	0	100
Creatinine	100	0	100

DOA Array Oral Fluid II

Repeatability

Control material prepared at the cut-off and ±50% cut-off was assessed. Each sample was assessed against the cut-off material twice a day for 10 days (n=80 for each sample). The percentage agreement was calculated for the number of samples that reported negative and positive correctly.

Assay	-50% cut-off	Cut-off	+50% cut-off	Agreement (%)
α-PVP	0	34	79	99.4
Amphetamine	0	36	80	100
Barbiturates	0	45	80	100
Benzodiazepines I	0	42	80	100
Benzodiazepines II	0	30	80	100
Benzoyllecgonine/Cocaine	0	38	80	100
Buprenorphine	0	37	79	99.4
Fentanyl	0	69	80	100
JWH-018	0	47	80	100
Ketamine	0	58	79	99.4
LSD	0	37	80	100
6-MAM	1	61	79	98.6
Metadone	0	74	80	100
Methamphetamine	0	59	80	100
Opiate	0	60	80	100
Oxycodone	1	15	80	99.4
PCP	0	37	79	99.4
Tramadol	0	73	80	100
THC	0	53	80	100
UR-144	0	12	80	100

Accuracy

Spiked samples at varying concentrations were assessed (50 spiked positive samples prepared at concentrations >cut-off, 10 negative spiked samples prepared at concentrations <cut-off and 40 blank negative samples). Each sample was assessed against the cut-off material to determine a positive or negative result. The percentage agreement was calculated as the percentage of correct reports out of the total number of analysed samples (n=100).

Assay	Spike +	Spike -	Agreement (%)
α-PVP	50	0	100
Amphetamine	50	0	100
Barbiturates	50	0	100
Benzodiazepines I	50	0	100
Benzodiazepines II	49	0	99
Benzoyllecgonine/Cocaine	49	0	99
Buprenorphine	49	0	99
Fentanyl	50	0	100
JWH-018	50	10	90
Ketamine	45	0	95
LSD	50	0	100
6-MAM	50	1	99
Metadone	50	0	100
Methamphetamine	50	0	100
Opiate	50	0	100
Oxycodone	49	0	99
PCP	49	0	99
Tramadol	50	0	100
THC	50	2	98
UR-144	50	0	100

Conclusion

Data indicate optimal analytical performance and applicability of the Evidence MultiSTAT system to the simultaneous screening of drugs of abuse in <20 minutes in urine and oral fluid.